

SPACA9 Antibody



PACO38710

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

cytoplasmic microtubule.

Gene ID:

SPACA9

Uniprot

Q96E40

Synonyms:

Sperm acrosome-associated protein 9, SPACA9, C9orf9

Immunogen:

Recombinant Human Sperm acrosome-associated protein 9 protein (1-222AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A