

XAGE2 Antibody



PACO38726

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

XAGE2

Uniprot

Q96GT9

Synonyms:

X antigen family member 2 (XAGE-2) (Cancer/testis antigen 12.2) (CT12.2) (G antigen family D member 3), XAGE2, GAGED3 XAGE2B

Immunogen:

Recombinant Human X antigen family member 2 protein (1-111AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A