

COG3 Antibody, Biotin conjugated



PACO38741

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Involved in ER-Golgi transport.

Gene ID:

COG3

Uniprot

Q96JB2

Synonyms:

Conserved oligomeric Golgi complex subunit 3 (COG complex subunit 3) (Component of oligomeric Golgi complex 3) (Vesicle-docking protein SEC34 homolog) (p94), COG3, SEC34

Immunogen:

Recombinant Human Conserved oligomeric Golgi complex subunit 3 protein (101-400AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A