

TBX15 Antibody, Biotin conjugated



PACO38769

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Probable transcriptional regulator involved in the development of the skeleton of the limb, vertebral column and head. Acts by controlling the number of mesenchymal precursor cells and chondrocytes.

Gene ID:

TBX15

Uniprot

Q96SF7

Synonyms:

T-box transcription factor TBX15 (T-box protein 15) (T-box transcription factor TBX14) (T-box protein 14), TBX15, TBX14

Immunogen:

Recombinant Human T-box transcription factor TBX15 protein (1-496AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A