

PRSS35 Antibody, HRP conjugated



PACO38787

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

PRSS35

Uniprot

Q8N3Z0

Synonyms:

Inactive serine protease 35, PRSS35, C6orf158

Immunogen:

Recombinant Human Inactive serine protease 35 protein (201-413AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A