

PRSS35 Antibody, FITC conjugated



PACO38788

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	PRSS35
Source:	Uniprot
Rabbit	Q8N3Z0
Isotype:	Synonyms:
IgG	Inactive serine protease 35, PRSS35, C6orf158
Applications:	Immunogen:
ELISA	Recombinant Human Inactive serine protease 35 protein (201-413AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A