

ADPRHL1 Antibody, HRP conjugated



PACO38815

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

ADPRHL1

Uniprot

Q8NDY3

Synonyms:

[Protein ADP-ribosylarginine] hydrolase-like protein 1 (EC 3.2. -. -) (ADP-ribosylhydrolase 2), ADPRHL1, ARH2

Immunogen:

Recombinant Human [Protein ADP-ribosylarginine] hydrolase-like protein 1 protein (1-272AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A