

# ADPRHL1 Antibody, FITC conjugated



PACO38816

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ADPRHL1

**Uniprot**

Q8NDY3

**Synonyms:**

[Protein ADP-ribosylarginine] hydrolase-like protein 1 (EC 3.2. -. -) (ADP-ribosylhydrolase 2), ADPRHL1, ARH2

**Immunogen:**

Recombinant Human [Protein ADP-ribosylarginine] hydrolase-like protein 1 protein (1-272AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A