

FAM45A Antibody



PACO38834

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

FAM45A

Uniprot

Q8TCE6

Synonyms:

Protein FAM45A, FAM45A

Immunogen:

Recombinant Human Protein FAM45A protein (201-357AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A