

NUDCD2 Antibody



PACO38866

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May regulate the LIS1/dynein pathway by stabilizing LIS1 with Hsp90 chaperone.

Gene ID:

NUDCD2

Uniprot

Q8WVJ2

Synonyms:

NudC domain-containing protein 2, NUDCD2

Immunogen:

Recombinant Human NudC domain-containing protein 2 protein (1-157AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A