

C1orf54 Antibody



PACO38870

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

C1orf54

Uniprot

Q8WWF1

Synonyms:

Uncharacterized protein C1orf54, C1orf54


Immunogen:

Recombinant Human Uncharacterized protein C1orf54 protein (17-131AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO38870 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO38870 at dilution of 1:100.