

PPDPFL Antibody



PACO38874

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

cell fate specification, exocrine pancreas development.

Gene ID:

PPDPFL

Uniprot

Q8WWR9

Synonyms:

Pancreatic progenitor cell differentiation and proliferation factor-like protein (Exocrine differentiation and proliferation factor-like protein), PPDPFL, C8orf22


Immunogen:

Recombinant Human Pancreatic progenitor cell differentiation and proliferation factor-like protein (1-81AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human ovarian cancer using PACO38874 at dilution of 1:100.