

CLEC4D Antibody, Biotin conjugated



PACO38885

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Functions as an endocytic receptor. May be involved in antigen uptake at the site of infection, either for clearance of the antigen, or for processing and further presentation to T cells.

Gene ID:

CLEC4D

Uniprot

Q8WXI8

Synonyms:

C-type lectin domain family 4 member D (C-type lectin superfamily member 8) (C-type lectin-like receptor 6) (CLEC-6), CLEC4D, CLECSF8 MCL

Immunogen:

Recombinant Human C-type lectin domain family 4 member D protein (39-215AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A