

# ANKRD7 Antibody



PACO38894

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

male gonad development.

**Gene ID:**

ANKRD7

**Uniprot**

Q92527

**Synonyms:**

Ankyrin repeat domain-containing protein 7 (Testis-specific protein TSA806), ANKRD7

**Immunogen:**

Recombinant Human Ankyrin repeat domain-containing protein 7 protein (54-254AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A