

PACO38918

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Co-chaperone that facilitates HSP-mediated activation of TSSK6.

Gene ID:

TSACC

Uniprot

Q96A04

Synonyms:

TSSK6-activating co-chaperone protein (SSTK-interacting protein) (SSTK-IP), TSACC, C1orf182


Immunogen:

Recombinant Human TSSK6-activating co-chaperone protein (1-125AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human gastric cancer using PACO38918 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human prostate cancer using PACO38918 at dilution of 1:100.