

RAB37 Antibody, FITC conjugated



PACO38924

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

endoplasmic reticulum-Golgi intermediate compartment, endosome, plasma membrane, secretory granule membrane, GTPase activity, intracellular protein transport, protein localization to plasma membrane, protein secretion, Rab protein signal transduction, regulation of exocytosis.

Gene ID:

RAB37

Uniprot

Q96AX2

Synonyms:

Ras-related protein Rab-37, RAB37

Immunogen:

Recombinant Human Ras-related protein Rab-37 protein (2-220AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A