

BORCS7 Antibody



PACO38930

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

Gene ID:

BORCS7

Uniprot

Q96B45

Synonyms:

BLOC-1-related complex subunit 7 (Diaskedin), BORCS7, C10orf32


Immunogen:

Recombinant Human BLOC-1-related complex subunit 7 protein (2-105AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO38930 at dilution of 1:100.



Immunofluorescent analysis of HeLa cells using PACO38930 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).