

ARL8A Antibody



PACO38934

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May play a role in lysosomes motility. Alternatively, may play a role in chromosome segregation.

Gene ID:

ARL8A

Uniprot

Q96BM9

Synonyms:

ADP-ribosylation factor-like protein 8A (ADP-ribosylation factor-like protein 10B) (Novel small G protein indispensable for equal chromosome segregation 2), ARL8A, ARL10B GIE2

Immunogen:

Recombinant Human ADP-ribosylation factor-like protein 8A protein (20-186AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A