

# CMTM5 Antibody, FITC conjugated



PACO38964

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:****Gene ID:**

CMTM5

**Uniprot**

Q96DZ9

**Synonyms:**

CKLF-like MARVEL transmembrane domain-containing protein 5 (Chemokine-like factor superfamily member 5), CMTM5, CKLFSF5

**Immunogen:**

Recombinant Human CKLF-like MARVEL transmembrane domain-containing protein 5 protein (1-35AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A