

CMTM5 Antibody, Biotin conjugated



PACO38965

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

CMTM5

Uniprot

Q96DZ9

Synonyms:

CKLF-like MARVEL transmembrane domain-containing protein 5 (Chemokine-like factor superfamily member 5), CMTM5, CKLFSF5

Immunogen:

Recombinant Human CKLF-like MARVEL transmembrane domain-containing protein 5 protein (1-35AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A