

# SPG11 Antibody, FITC conjugated



PACO38992

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

cytoplasm, cytoplasmic vesicle, extracellular exosome, lysosomal membrane, nucleolus, plasma membrane.

**Gene ID:**

SPG11

**Uniprot**

Q96J17

**Synonyms:**

Spatacsin (Colorectal carcinoma-associated protein) (Spastic paraplegia 11 protein), SPG11, KIAA1840

**Immunogen:**

Recombinant Human Spatacsin protein (2019-2318AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A