

C9orf72 Antibody, HRP conjugated



PACO38999

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May play a role in endosomal trafficking and autophagy.

Gene ID:

C9orf72

Uniprot

Q96LT7

Synonyms:

Guanine nucleotide exchange C9orf72, C9orf72

Immunogen:

Recombinant Human Guanine nucleotide exchange C9orf72 protein (138-430AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A