

# C9orf72 Antibody, FITC conjugated



PACO39000

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May play a role in endosomal trafficking and autophagy.

**Gene ID:**

C9orf72

**Uniprot**

Q96LT7

**Synonyms:**

Guanine nucleotide exchange C9orf72, C9orf72

**Immunogen:**

Recombinant Human Guanine nucleotide exchange C9orf72 protein (138-430AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A