

TCEANC2 Antibody



PACO39014

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:**Gene ID:**

TCEANC2

Uniprot

Q96MN5

Synonyms:

Transcription elongation factor A N-terminal and central domain-containing protein 2, TCEANC2, C1orf83


Immunogen:

Recombinant Human Transcription elongation factor A N-terminal and central domain-containing protein 2 protein (1-208AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human bladder cancer using PACO39014 at dilution of 1:100.