

GLT1D1 Antibody



PACO39018

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

cytoplasm.

Gene ID:

GLT1D1

Uniprot

Q96MS3

Synonyms:

Glycosyltransferase 1 domain-containing protein 1 (EC 2.4. -. -), GLT1D1


Immunogen:

Recombinant Human Glycosyltransferase 1 domain-containing protein 1 protein (17-152AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human bladder cancer using PACO39018 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human thyroid tissue using PACO39018 at dilution of 1:100.