

# PRELID3A Antibody



PACO39026

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

mitochondrial intermembrane space, phosphatidic acid, transporter activity, phospholipid transport, positive regulation of protein targeting to mitochondrion.

**Gene ID:**

PRELID3A

**Uniprot**

Q96N28

**Synonyms:**

PRELI domain containing protein 3A (Protein slowmo homolog 1), PRELID3A, C18orf43 SLMO1

**Immunogen:**

Recombinant Human PRELI domain containing protein 3A protein (1-172AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A