

RASL10B Antibody, Biotin conjugated



PACO39053

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May facilitate the release of atrial natriuretic peptide by cardiomyocytes and hence play a role in the regulation of arterial pressure.

Gene ID:

RASL10B

Uniprot

Q96S79

Synonyms:

Ras-like protein family member 10B (Ras-like protein VTS58635) (Ras-related protein 17) (RRP17), RASL10B

Immunogen:

Recombinant Human Ras-like protein family member 10B protein (15-200AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A