

CNPY4 Antibody



PACO39058

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Plays a role in the regulation of the cell surface expression of TLR4.

Gene ID:

CNPY4

Uniprot

Q8N129

Synonyms:

Protein canopy homolog 4, CNPY4


Immunogen:

Recombinant Human Protein canopy homolog 4 protein (22-248AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO39058 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human appendix tissue using PACO39058 at dilution of 1:100.