

ZFP82 Antibody, HRP conjugated



PACO39063

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZFP82

Uniprot

Q8N141

Synonyms:

Zinc finger protein 82 homolog (Zfp-82) (Zinc finger protein 545), ZFP82, KIAA1948
ZNF545

Immunogen:

Recombinant Human Zinc finger protein 82 homolog protein (1-200AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A