

# ZFP82 Antibody, Biotin conjugated



PACO39065

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May be involved in transcriptional regulation.

**Gene ID:**

ZFP82

**Uniprot**

Q8N141

**Synonyms:**

Zinc finger protein 82 homolog (Zfp-82) (Zinc finger protein 545), ZFP82, KIAA1948  
ZNF545

**Immunogen:**

Recombinant Human Zinc finger protein 82 homolog protein (1-200AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A