

D2HGDH Antibody, HRP conjugated



PACO39079

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalyzes the oxidation of D-2-hydroxyglutarate to alpha-ketoglutarate.

Gene ID:

D2HGDH

Uniprot

Q8N465

Synonyms:

D-2-hydroxyglutarate dehydrogenase, mitochondrial (EC 1.1.99. -), D2HGDH, D2HGD

Immunogen:

Recombinant Human D-2-hydroxyglutarate dehydrogenase, mitochondrial protein (14-521AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A