

SAMD3 Antibody



PACO39094

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

SAMD3

Uniprot

Q8N6K7

Synonyms:

Sterile alpha motif domain-containing protein 3 (SAM domain-containing protein 3), SAMD3


Immunogen:

Recombinant Human Sterile α motif domain-containing protein 3 protein (62-197AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human pancreatic tissue using PACO39094 at dilution of 1:100.