

PNPLA1 Antibody, HRP conjugated



PACO39107

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Lipid hydrolase. Important in the formation of the epidermal lipid barrier. Plays a role in glycerophospholipid metabolism.

Gene ID:

PNPLA1

Uniprot

Q8N8W4

Synonyms:

Patatin-like phospholipase domain-containing protein 1 (EC 3.1.1. -), PNPLA1

Immunogen:

Recombinant Human Patatin-like phospholipase domain-containing protein 1 protein (101-300AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A