

# FSIP1 Antibody, HRP conjugated



PACO39119

---

## Product Information

|                               |  |
|-------------------------------|--|
| <b>Size:</b>                  | <b>Protein Background:</b>   |
| 50ug                          |  |
| <b>Reactivity:</b>            | <b>Gene ID:</b>  |
| Human                         | FSIP1  |
| <b>Source:</b>                | <b>Uniprot</b>   |
| Rabbit                        | Q8NA03   |
| <b>Isotype:</b>               | <b>Synonyms:</b>   |
| IgG                           | Fibrous sheath-interacting protein 1, FSIP1                                    |
| <b>Applications:</b>          | <b>Immunogen:</b>  |
| ELISA                         | Recombinant Human Fibrous sheath-interacting protein 1 protein (101-400AA).    |
| <b>Recommended dilutions:</b> | <b>Storage:</b>  |
|                               | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |

## Product Images

---

N/A

N/A