

FSIP1 Antibody, Biotin conjugated



PACO39121

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

FSIP1

Uniprot

Q8NA03

Synonyms:

Fibrous sheath-interacting protein 1, FSIP1

Immunogen:

Recombinant Human Fibrous sheath-interacting protein 1 protein (101-400AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A