

PACO39122

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000

Protein Background:

cytoplasmic membrane-bounded vesicle.

Gene ID:

SPERT

Uniprot

Q8NA61

Synonyms:

Spermatid-associated protein (Protein chibby homolog 2), SPERT, CBY2


Immunogen:

Recombinant Human Spermatid-associated protein (1-412AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: SPERT antibody at 7 μ g/ml. Lane 1: Rat gonadal tissue. Lane 2: Mouse lung tissue. Lane 3: A549 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 52, 48 kDa. Observed band size: 52 kDa.