

ZNF738 Antibody



PACO39142

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZNF738

Uniprot

Q8NE65

Synonyms:

Protein ZNF738, ZNF738


Immunogen:

Recombinant Human Protein ZNF738 protein (1-137AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human liver tissue using PACO39142 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human lung tissue using PACO39142 at dilution of 1:100.