

GPR156 Antibody



PACO39158

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Orphan receptor.

Gene ID:

GPR156

Uniprot

Q8NFN8

Synonyms:

Probable G-protein coupled receptor 156 (G-protein coupled receptor PGR28) (GABAB-related G-protein coupled receptor), GPR156, GABABL PGR28


Immunogen:

Recombinant Human Probable G-protein coupled receptor 156 protein (614-814AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO39158 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human brain tissue using PACO39158 at dilution of 1:100.