

NFKBID Antibody, HRP conjugated



PACO39171

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May regulate the expression of IL-2, IL-6, and other cytokines through regulation on NF-kappa-B activity. Functions in the regulation of inflammatory responses. May also regulate TCR-induced negative selection of thymocytes.

Gene ID:

NFKBID

Uniprot

Q8NI38

Synonyms:

NF-kappa-B inhibitor delta (I-kappa-B-delta) (IkB-delta) (IkappaBdelta) (IkappaBNS) (T-cell activation NFKB-like protein) (TA-NFKBH), NFKBID, IKBNS

Immunogen:

Recombinant Human NF-kappa-B inhibitor delta protein (1-150AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A