

# NFKBID Antibody, FITC conjugated



PACO39172

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May regulate the expression of IL-2, IL-6, and other cytokines through regulation on NF-kappa-B activity. Functions in the regulation of inflammatory responses. May also regulate TCR-induced negative selection of thymocytes.

**Gene ID:**

NFKBID

**Uniprot**

Q8NI38

**Synonyms:**

NF-kappa-B inhibitor delta (I-kappa-B-delta) (IkB-delta) (IkappaBdelta) (IkappaBNS) (T-cell activation NFKB-like protein) (TA-NFKBH), NFKBID, IKBNS

**Immunogen:**

Recombinant Human NF-kappa-B inhibitor delta protein (1-150AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A