

Tmprss4 Antibody, FITC conjugated



PACO39196

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Probable protease. Seems to be capable of activating ENaC.

Gene ID:

Tmprss4

Uniprot

Q8VCA5

Synonyms:

Transmembrane protease serine 4 (EC 3.4.21. -) (Channel-activating protease 2) (mCAP2), Tmprss4, Cap2

Immunogen:

Recombinant Mouse Transmembrane protease serine 4 protein (52-435AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A