

# UBTD2 Antibody, FITC conjugated



PACO39208

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

UBTD2

**Uniprot**

Q8WUN7

**Synonyms:**

Ubiquitin domain-containing protein 2 (Dendritic cell-derived ubiquitin-like protein) (DC-UbP) (Ubiquitin-like protein SB72), UBTD2, DCUBP

**Immunogen:**

Recombinant Human Ubiquitin domain-containing protein 2 protein (45-234AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A