

TCTEX1D2 Antibody



PACO39214

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

TCTEX1D2

Uniprot

Q8WW35

Synonyms:

Tctex1 domain-containing protein 2, TCTEX1D2


Immunogen:

Recombinant Human Tctex1 domain-containing protein 2 protein (1-142AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO39214 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human ovarian cancer using PACO39214 at dilution of 1:100.