

SSMEM1 Antibody, HRP conjugated



PACO39219

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

SSMEM1

Uniprot

Q8WWF3

Synonyms:

Serine-rich single-pass membrane protein 1, SSMEM1, C7orf45

Immunogen:

Recombinant Human Serine-rich single-pass membrane protein 1 protein (56-244AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A