

MRPS31 Antibody, Biotin conjugated



PACO39245

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

mitochondrial inner membrane, mitochondrion, nucleolus, poly(A) RNA binding, protein domain specific binding, mitochondrial translation, mitochondrial translational elongation, mitochondrial translational initiation, mitochondrial translational termination, organelle organization.

Gene ID:

MRPS31

Uniprot

Q92665

Synonyms:

28S ribosomal protein S31, mitochondrial (MRP-S31) (S31mt) (Imogen 38)
(Mitochondrial small ribosomal subunit protein mS31), MRPS31, IMOGN38

Immunogen:

Recombinant Human 28S ribosomal protein S31, mitochondrial protein (66-395AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A