

DIRC1 Antibody




PACO39258

Product Information

Size: 50ug	Protein Background:
Reactivity: Human	Gene ID: DIRC1
Source: Rabbit	Uniprot Q969H9
Isotype: IgG	Synonyms: Disrupted in renal carcinoma protein 1 (Disrupted in renal cancer protein), DIRC1
Applications: ELISA, IHC	Immunogen: Recombinant Human Disrupted in renal carcinoma protein 1 protein (1-104AA).
Recommended dilutions: ELISA:1:2000-1:10000, IHC:1:20-1:200	Storage: Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO39258 at dilution of 1:100.