

PLEKHB2 Antibody, HRP conjugated



PACO39279

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in retrograde transport of recycling endosomes.

Gene ID:

PLEKHB2

Uniprot

Q96CS7

Synonyms:

Pleckstrin homology domain-containing family B member 2 (PH domain-containing family B member 2) (Evectin-2), PLEKHB2, EVT2

Immunogen:

Recombinant Human Pleckstrin homology domain-containing family B member 2 protein (1-222AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A