

PACO39310

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Co-chaperone that binds directly to HSC70 and HSP70 and regulates their ATPase activity.

**Gene ID:**

SGTB

**Uniprot**

Q96EQ0

**Synonyms:**

Small glutamine-rich tetratricopeptide repeat-containing protein beta (Beta-SGT) (Small glutamine-rich protein with tetratricopeptide repeats 2), SGTB, SGT2

**Immunogen:**


Recombinant Human Small glutamine-rich tetratricopeptide repeat-containing protein &beta; protein (14-304AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human colon tissue using PACO39310 at dilution of 1:100.