

SCMH1 Antibody



PACO39322

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Associates with Polycomb group (PcG) multiprotein complexes; the complex class is required to maintain the transcriptionally repressive state of some genes.

Gene ID:

SCMH1

Uniprot

Q96GD3

Synonyms:

Polycomb protein SCMH1 (Sex comb on midleg homolog 1), SCMH1

Immunogen:

Recombinant Human Polycomb protein SCMH1 protein (201-500AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A