

# ZC3H10 Antibody, Biotin conjugated



PACO39333

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

poly(A) RNA binding.

**Gene ID:**

ZC3H10

**Uniprot**

Q96K80

**Synonyms:**

Zinc finger CCCH domain-containing protein 10, ZC3H10, ZC3HDC10

**Immunogen:**

Recombinant Human Zinc finger CCCH domain-containing protein 10 protein (1-300AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A